



Experience ultra-smooth insertion with our Transitionless-Tip™.



10 cm Transitionless-Tip
introducer catheter



40 cm nitinol wire guide
with palladium tip



21 gage, 7 cm EchoTip®
needle

Micropuncture®

INTRODUCER SET

Fewer complications

In a study by Shah et al, access-site bleeding complication rates dropped from 11.7% to 2.6% when Micropuncture was used.

Smooth entry

The introducer's Transitionless-Tip™ has been engineered to glide into the vessel with less force.

High visibility

The sharp, 21 gage needle features EchoTip® technology that enhances ultrasound visibility. The supportive nitinol wire has a palladium tip for excellent visibility under fluoroscopy.



cookmedical.com